



Fog droplets can settle under static conditions and can keep floating in turbulent conditions



It is recommended to charge once every three months if it has not been used for a long time

# PHARMA 20

## LITHIUM STERILIZER SPRAYER

Mist sprayer, operated by lithium batteries, is ergonomically designed for long hours of comfortable use. This is very well suited for use in any situations that require rapid sanitizing such as the current Covid-19 pandemic. Sanitizing is made affordable with the Pharma 20 Lithium Sterilizer Sprayer.



**VOLTAGE**  
12 V



**CURRENT**  
2.0 A



**BATTERY**  
2000 Ma



**CAPACITY**  
5 L



**OPERATING HOURS**  
Approx 0.5 Hr



**UNLOADED WEIGHT**  
1.5 Kg



**LOADED WEIGHT**  
7.5 Kg



**CHARGING TIME**  
1 Hr



**BATTERY  
OPERATING TIME**  
1 Hr



**SPRAY DISTANCE**  
1.5-2.0 M



**DROPLET SIZE**  
≥ 50 Um



**ATOMIZATION FLOW**  
≥ 142-145 ml / Min

### SCOPE OF APPLICATION:

Airport, hotel, hospital, bar, bus, supermarket, movie theater, children's playground.

### PRECAUTIONS:

Ensure that the liquid is clean and free of impurities. If the amount of mist becomes smaller, you can unscrew the nozzle to clean it. After the container is filled, the machine will have 30 seconds of pipe exhaust time and then work normally.